

ROEQ

ROBOTIC EQUIPMENT

We move your business

ROEQ R150 Rack

The R150 is developed to support internal logistics in industry and health care facilities.

The unit is made as a base and can be used with a standard rack attached or as a base for a customized application.



Application interface

The 4 corner holes are prepared for attaching various end user applications.

Frame size

800 x 800 mm

Maximum payload

MiR100: 50 kg
MiR200: 150 kg

Top Module Rack (TMR)

The Top Module for Racks is made with a click-in and lifting mechanism that enables the system to load and unload Racks.

Docking

When MiR100/MiR200 places the Rack in the Docking Station, the Rack will be aligned to e.g. an assembly line.

This is how the ROEQ R150 works



MiR100/MiR200 with R150 approaches the Docking Station.

MiR100/MiR200 releases the Rack.

ROEQ R150 stays attached to the Docking Station – MiR100/MiR200 is off on a new mission.

Visit our website for more information

www.roeq.dk

VICOSYSTEMS